

Product Specifications

Brand : Klinion®
Product : Kliniderm® Film Roll

Product description

Kliniderm® Film Roll is made of polyurethane foil with a skin friendly acrylic adhesive. The films are semi-permeable to exchange oxygen and water vapor but to prevent microorganisms to enter the wound.

Intended use

Wound film roll for surface wounds with slight exudate.

Classification

Class I under the Medical Devices Regulations (SI 2002 No. 618, as amended) and
Class I under Regulation (EU) 2017/745 Medical Device Regulation (MDR).

Packaging & Product Reference

Product Reference	Description	Size (*)	Piece(s)	Inner Carton (pcs)	Outer Carton (pcs)
40514870	Kliniderm Film Roll 5 cm x 10 m	5 cm x 10 m	1	1	50
40514871	Kliniderm Film Roll 10 cm x 10 m	10 cm x 10 m	1	1	20
40514872	Kliniderm Film Roll 15 cm x 10 m	15 cm x 10 m	1	1	20

Product Reference	Description	Size (*)	Piece(s)	EAN Inner Carton	EAN Outer Carton
40514870	Kliniderm Film Roll 5 cm x 10 m	5 cm x 10 m	1	8715343015674	8715343015681
40514871	Kliniderm Film Roll 10 cm x 10 m	10 cm x 10 m	1	8715343015698	8715343015711
40514872	Kliniderm Film Roll 15 cm x 10 m	15 cm x 10 m	1	8715343015728	8715343015742

*Sizes are approximate

Features

- Transparent
- Semi-permeable
- Self-adhesive

Product Specifications

Precautions

- Single-use
- Non-sterile
- Store at room temperature, protected from direct sunlight.

Specifications

Component	Material	CAS number
Film	PU Film (17±5 g/m2)	78-93-3 7631-86-9 68-12-2
Glue	Acrylic adhesive (≥25 g/m2)	141-78-6 103-11-7
Protective layers	PET Film (0.025±0.005 mm)	63148-62-9 25038-59-9
	Glassine paper (62±5 g/m2)	65996-61-4 7732-18-5 2627-95-4

REF number	Size	Holding power	180° peel strength	Water vapor transmission rate
40514870	5±0.5cm x 10±0.5m	≤2.5 mm	≥1N/cm	≥1000g/m2/24h
40514871	10±0.5cm x 10±0.5m	≤2.5 mm	≥1N/cm	≥1000g/m2/24h
40514872	15±0.5cm x 10±0.5m	≤2.5 mm	≥1N/cm	≥1000g/m2/24h

Attributes	Test method/Standard
Holding power	GB/T4851-2014 - Test Method for Holding Power of Pressure Sensitive Adhesive Tapes
180° Peel strength	GB/T2792-2014 - Pressure Sensitive adhesive tape peel strength test method
Water vapor transmission rate	YY/T0471.2-2004 - Water Vapor Transmission Rate of Transparent Dressing Film

Product Specifications

Presence of specific components

Component	Present
Latex	No
Phthalates	No
Substances of animal origin	No
Substances of human origin	No
Blood or blood derivates	No
Medicinal substances	No
Heavy metals	No
PVC	No

Sterilization

Non-sterile.

Storage conditions

Kliniderm® Film Roll should be stored at room temperature, protected from direct sunlight.

Shelf life

Kliniderm® Film Roll has a shelf life of 5 years from production date, when stored under specified conditions.

Picture

Product Specifications



Product ref.: 40514871

Approval	Name	Date	Signature
Product Manager Wound Care	Florin Lazar	23.05.2025	
RA Specialist	Emil Laske	23.05.2025	